

How to copy recipes on the EVG 810 Plasma Activation

1. Select 'Login Here' (lower left corner of the screen) and use the pop-up window to login to the system.
2. Verify you are in the correct recipe directory by looking at the banner at the top of the program window, it with read EVGVISU – *directory name*.ecp. The general use directory is called 'Official Recipe Directory.ecp'. If you need to open a different directory, perform the following:
 - a. Select 'Recipe'.
 - b. Select File → Open. A pop-up window will appear that should be in the 'Recipes' folder, double click on a directory to load it.

NOTE: If you are not in the recipes folder you will need to locate it, use the drop down menu in the pop-up window to begin looking for it. The file path is:
C:\Program Files\EVG\EVGVisu\Recipes
3. Left click on the recipe you want to copy to highlight it.
4. Right click and select 'Copy Recipe' from the pop-up window.
5. If you are copying the recipe to the same directory:
 - a. Right click again and select 'Paste Recipe' from the pop-up window. The copied recipe will appear in the directory with the name 'Copy of *recipe name*'.
6. If you are copying the recipe to another directory:
 - a. Open the directory per Step 2.b.
 - b. Left click on any recipe in the directory to highlight it.
 - c. Right click and select 'Paste Recipe' from the pop-up window. The copied recipe will appear in the directory with the name 'Copy of *recipe name*'.
7. To edit the name of the copied recipe:
 - a. Left click on the recipe to highlight it.
 - b. Left click again to enable name editing.
 - c. Use the keyboard to edit the name and press 'Enter' when finished.
8. Left click on the save icon (floppy disk) at the top of the screen to save your changes.